

# How to Support a 30 Degree Turned Position Using The Backbone™ Positioning Wedge

The BackBone is quick and easy to place

1. Roll patient as far as possible to side position



2. Angle The BackBone so that “tips” point toward bed



3. Press upper tip under scapula



4. Press Lower tip under buttocks



5. Place to avoid sacrum, or wound area



6. Have your patient roll back onto The BackBone



Place directly under patient or under moisture control pads/chux. Always place ABOVE / ON TOP OF patient positioning sheets (above sliders or draw sheets)